
ROXX

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## Contents

*This manual has been prepared with the utmost care and is based on information known to Life & Mobility at the time of publication.*

*The information in this manual applies to the standard version of the product. Life & Mobility does not accept responsibility for any damage that may occur due to changes or modifications to the product after selling it for the first time. Neither can Life & Mobility be held liable for any damage caused by third-party work.*

*The information in this manual may not be reproduced for commercial goals and/or made public by any means without prior written permission from Life & Mobility.*

*Pursuant to the Trade Names Act, the trade name and trademarks used by Life & Mobility may not be considered free.*

*No rights can be derived from this text. Subject to changes.*

*© Life & Mobility – All rights reserved*

Foreword	4
Product description	5
Matrixx seat system	6
Options	7
Technical specifications	8
Safety instructions	9
Settings and adjustments	14
Folding the wheelchair	26
Use	28
Transport	30
Maintenance	33
Technical problems	34
Guarantee	35
Used wheelchairs/environment	36

## Foreword

The Roxx has a modular, multi-adjustable frame and can accommodate many sizes, making it possible to accurately adjust it to the individual. In combination with the correct seat support from the Matrixx seat system, you can create a perfectly fitting wheelchair for every user.

The Roxx is splash-resistant and is suitable for both indoor and outdoor use. The maximum occupant weight is 160 kg. This is indicated on the CE sticker.

The Roxx has passed the tests confirming compliance with the relevant requirements of the European wheelchair standards NEN-EN 12182 and NEN-EN 12183. It has also shown to be in compliance with the associated crashtest safe wheelchair standards in ISO 7176-19:2008 (see page 30). The Roxx therefore complies with the mandatorily applicable European REGULATION (EU) 2017/745 for Medical Devices.

**!** **Important:** before using the Roxx, please read these instructions carefully.

**!** **Important:** failing to follow the safety regulations, instructions, warnings and maintenance recommendations, can negatively affect the operation of the Roxx and possibly result in physical injury or damage to the wheelchair or the environment.

Sections in this document that are relevant to your safety are marked with a safety symbol.

## Product description

### Delivery Check

Before you start using it, check that the wheelchair delivered matches the product description below. If you notice any damage caused during transport, please inform your dealer. We at Life & Mobility are sure that your Roxx will meet your expectations and we wish you much success with it.

### Product description

1. Push-handle
2. Matrixx backrest
3. Side with armrest
4. Matrixx seat
5. Matrixx legrest
6. Brake
7. Folding or rigid frame
8. Front wheel
9. Back wheel
10. Headrest adapter
11. Anti-tipping support (optional)
12. Front and back taxi hooks on both left and right

### Not illustrated

- User manual
- Options/ Accessories

These instructions describe the standard products. The Roxx can of course be extended with many options and accessories, these are listed under options.

## Matrixx seat system

Optimal seating support is achieved by a specific combination of a basic wheelchair and the correct seat elements. The seat and backrest angle of a Life & Mobility wheelchair can be perfectly adjusted to suit any user. In addition, the various seat elements can be chosen modularly within the Matrixx seat system, dependent on the functional limitations of the user. This allows a suitable wheelchair to be created for any user.

### Advantages of the Matrixx seat system

- Interchangeability of elements from the Life & Mobility wheelchairs
- Modular product structure with many options

### Advantages for the user

- Optimal seating support
- Good pressure distribution
- Perfect sitting position
- Excellent sitting comfort

### The Matrixx seat system consists of the following elements

- Seats
- Backrests
- Legrests
- Armrests
- Headrests
- Seat elements
- Matrixx cleaner

For more information and the latest news about the Matrixx seat system, please consult the Life & Mobility website [www.life-mobility.com](http://www.life-mobility.com)

## Options

In addition to the standard components and the Matrixx seat system, the following options are available for the Roxx:

### Frame

- Folding or rigid frame

### Swivel castors

- Semi-solid or solid tyres

### Rear wheels

- Semi-solid tyres or pneumatic tyres
- Single-hand drive
- Hand rim covers
- Transparent spoke protectors

### Brake system

- Standard or extended brake
- Drum brakes operated by the accompanying person
- Single-hand brake system
- Hub brakes with brake levers

### Wheel plates

- Standard or extended wheel plates

### Backrest angle

- Once-only or continuous adjustment

### Sides

- Removable or lifting sides
- Extended side shields

### Accessories

- Worktop
- Anti-tipping support
- Stick holder
- Seat belt

For more information and the latest news about the various options, please consult the Life & Mobility website [www.life-mobility.com](http://www.life-mobility.com)

## Technical specifications

Seat width	37 to 60 cm (in steps of 3 cm)
Seat depth	42 to 58 cm
Seat height	38 to 60 cm (in steps of 1 cm)
Armrest height	18 to 27 cm (in steps of 1,5 cm)
Width of wheelchair	60 to 80 cm (in steps of 5 cm)
Width of wheelchair folded (incl. wheels)	31 cm

Seat angle	0° or adjustable
Backrest angle	Once-only adjustment (0 – 180°) Continuously adjustable (85 – 105° / 115°)
Backrest height	Adjustable in 1 cm steps
Camber angle	2°
Lower leg length	35 - 45 cm / 40 - 55 cm / 50 - 65 cm
Total weight	From 25 kg
Lift weight	From 10 kg
Max. weight of occupant	160 kg

Complies with EN 12182 / EN 12183 / ISO 7176-19:2008

Life & Mobility recommends using the Roxx in an ambient temperature of between -5 °C and +40 °C.

### Service and technical support

Consult your dealer if you have questions about specific adjustments, maintenance or repair. Always have the serial number available. You can find this on the CE sticker on the wheelchair (A).

### Expected service life

Life & Mobility assigns a service life of 7 years to this product. This means that during this period, the product is fixable, safe and usable.

The statement above is only applicable if the user uses the product in the manner that Life & Mobility intended, marked with CE and including any original accessories.

A

## Safety instructions

Life & Mobility accepts no liability for damage or injury caused by negligence or non-compliance with the safety instructions.

Dependent on the specific circumstances or accessories used, supplementary safety instructions can be required.

Allways notify Life & Mobility and the relevant authority in case of a serious incident.

### General

- Never change the construction of the wheelchair.
- Only use original Life & Mobility parts and accessories. These are tested and safe.

### Markings on the wheelchair

- Markings, safety symbols and instructions are displayed on this wheelchair. They must not be covered or removed. They must remain on and be clearly visible throughout the life of the wheelchair.

### Assembly and adjustment

- Before using the wheelchair, check that all of the components of the wheelchair are correctly and securely installed.
- Ensure that the push-handles, the backrest height, the armrests, the lower leg length and the headrest, if present, are adjusted to the right height.
- Ensure that there is a minimum distance of 5 cm between the legrests and the floor.
- Check the dimensions and correctly adjust them in the manner described in this manual.

### Sitting in the wheelchair

- Ensure that the fingers and limbs of the occupant do not become trapped. Be particularly vigilant when:
  - The occupant is sitting lopsidedly in the wheelchair.

## Safety instructions

- The armrests are lifted backwards.
- The occupant does not have complete control of his/her arms and/or legs.

### Using the wheelchair

(see also page 28 for more information)

 **Note:** The Roxx is not fitted as standard with anti-tipping support. Check whether the wheelchair can tip over backwards if tilted suddenly backwards.

 **Note:** When travelling up a slope, the centre of gravity will move backwards, making it easier for the wheelchair to tip over backwards.

 **Note:** When travelling down a slope, the centre of gravity will move forwards, making it easier for the occupant to fall out of the wheelchair.

 **Note:** As a wheelchair user you are vulnerable in traffic. Ensure that you are visible to other road users. Observe the applicable traffic regulations.

 **Note:** Avoid infrequently used routes, ensure that in case of emergency help can be provided.

- In the following cases, the occupant or accompanying person should be alert to the possibility of the wheelchair tipping over backwards:
  - When the armrests are lifted backwards.
  - When the wheel balance/ ride adjustment is active.
  - When a bag is suspended from the push-handle.
  - When the backrest angle is large.
  - When travelling up a slope.
  - When someone is standing on the tipping lever and/ or the anti-tipper.

## Safety instructions

- When self propelling the wheelchair, consider your hands if there is not much space at the side(s) of the wheelchair.
- Move forward by means of the hoop at all times. If the user cannot do this for specific reasons, Life & Mobility recommends dismantling the brake per side and placing the brake assembly profile with the furthest protruding part back to the front. In doing so, assess whether the braking effect after conversion is sufficient and whether it does not conflict with any other options.
- Ensure that the removable sides are folded downwards before you start to move.
- Accompanying persons must never leave the wheelchair on a slope without engaging the brake.
- Never tackle obstacles higher than 10 cm on your own. Only experienced users with sufficient strength are generally capable of tackling obstacles up to 10 cm in height on their own.
- Engage the brake on the wheelchair before getting out of the chair or transferring from it.
- Only use the brake to lock the stationary wheels. The brake is not intended for slowing the wheelchair down.
- When braking, heat is generated by the friction between your hands and the hand rims. If you brake too strongly, your hands can become too hot.

### Transporting the wheelchair without the occupant

(see also page 26 and 27 for more information)

The Roxx equipped with a folding frame is easy to transport in a car. When doing so, consider the following points:

- Ensure that your fingers do not get caught when folding up the wheelchair.
- Always carry the wheelchair by the frame, never by the legrest or armrests.
- Ensure that the wheelchair is completely unfolded before you use it.

## Safety instructions

### Transporting the wheelchair with the occupant in it in a taxibus (see also page 30 for more information)

The Roxx was designed to be secured in a forward direction in a taxi van. The Roxx was tested positively as a crash-test safe wheelchair in accordance with ISO 7176-19:2008. It can therefore be transported with the occupant in it in a taxibus (see also the sticker on the wheelchair).

- Attach the wheelchair with its occupant in it in a taxibus with a four-point securing system.
- In addition, use a 3-point safety belt to safely transport the occupant of the wheelchair.

**Note:** The normal seat belt used to hold the occupant in position is not intended for use as a safety belt.

When transporting the wheelchair with the occupant in it, Life & Mobility recommends fitting the wheelchair with a Matrixx headrest.

### Slopes and other obstacles

(see also page 29 for more information)

- The wheelchair is more likely to tip over on a slope than on a level surface.
- Always go off a pavement backwards to prevent the wheelchair from leaning forwards allowing the occupant to fall out. Never go down stairs in the wheelchair alone.
- Never propel the wheelchair up a slope with an angle greater than 5°.
- Ask somebody to help you if you are unsure whether you can reach the top of the slope by yourself.
- Lean or push your body in the direction of the top of the slope. This will shift the centre of gravity of the wheelchair forward and reduce the risk of the chair tipping over backwards.

## Safety instructions

- The person pushing the wheelchair must continue to push the chair to prevent it from starting to move backwards under its own momentum.
- Do not turn the wheelchair round on a steep slope. If the wheelchair is facing across a steep slope, it can very easily tip over.
- Keep your speed under control when descending slopes.
- The accompanying person must reduce the speed of the wheelchair using the push-handles.
- If you are using the wheelchair yourself, slow the wheelchair down using your hands on the hand rims only. Never use the brake as this can cause you to lose control of the wheelchair.

### Other

- Take care when using the wheelchair on surfaces that are slippery due to rain, ice or snow!
- Do not allow the wheelchair to come in contact with seawater: seawater is aggressive and will corrode it.
- Do not allow the wheelchair to come in contact with sand: sand can penetrate the rotating parts of the wheelchair, causing unnecessarily rapid wear.
- Do not leave the wheelchair in direct sunlight: the upholstery can become too hot, causing burns or oversensitivity of the skin.
- If you are under the influence of substances that can influence your ability to drive, you are not allowed to use the wheelchair.
- Your eyesight must be sufficiently good to allow you to use the wheelchair safely in the situation concerned.
- Ensure that there is no loose-hanging clothing. It could become caught between the wheels.

## Settings and adjustments

### Settings and adjustments

The Roxx has a wide range of settings and adjustment possibilities that can be used to adjust the wheelchair to meet all of the requirements of the user. Life & Mobility recommends that you consult your dealer with regard to the correct settings. You can easily make adjustments yourself, without the use of tools.

#### Frame

##### Seat width

The width of the seat frame is adjustable via the frame and can be set to widths of 37, 42, 47, 52 and 57 cm. In addition, it is also possible to extend the sides that are fitted with armrests by 1.5 cm on each side. When changing the frame width, you need to install a different seat and backrest.

##### Adjusting the seat width via the frame




1. Loosen the 6 socket screws (A).
2. Set the desired seat width, ensure that the seat is widened symmetrical on both sides.
3. Retighten the 6 socket screws (A) (with the nut on the bottom).

##### Extending the sides outwards




1. Lift the operating lever (B1) and remove the side.

or;

1. Unscrew the socket screw (B2) of the lifting armrest a couple of turns and remove the side.

## Settings and adjustments

2. Using 2 Allen keys, unscrew the book screws (C) of the side guard.
3. Unscrew the 2 socket screws at the bottom of the side tube (D).
4. Turn over the plate (E).
5. Position the oval tube outwards on the plate.
6. Retighten the 2 socket screws at the bottom of the tube (D).
7. Retighten the book screw (C) of the side guard on the other side.
8. Reposition the entire side on the side frame.
9. Repeat this on the other side, if required.

**Note:** The side is locked when the operating lever has audibly clicked into the locked position and the side can no longer be removed.

#### Fitting the side pieces

1. Insert a M6x8 socket head screw in the threaded hole (F) and tighten firmly using an Allen key.
2. The side piece is now fixed securely and can no longer be folded down.
3. Repeat this process with the other side if necessary.

#### Seat depth

The seat depth can be set by replacing the 4 connecting pieces in the side frames.


1. Remove the sides.
2. Remove, if necessary, the seat and the brake.
3. Remove the connecting pieces by unscrewing the socket screws (A).
4. Position the desired connecting pieces 0-4-8-12 cm in the side frame.
5. Retighten the socket screws (A).

## Settings and adjustments

The seat depth can be set more accurately by shifting the backrest tubes by 4 cm.

1. Remove the seat.
2. Unscrew the 2 socket screws on each side (A) a few turns. We recommend you use the socket wrench M5 for this, as pictured.
3. Slide the plate with the push-handles forward or backward until the desired seat depth has been reached.
4. Determine the precise size using the graduation marks on both sides of the seat frame. These are indicated using stickers.
5. If necessary, adjust the ride characteristics.

**Note:** Ensure that the depth is the same at the left and right.


6. Retighten the socket screws (A).

**Note:** Ensure that the side socket (B) connects correctly to the side frame.

### Adjusting the axle plate

1. Remove the wheel.
2. Unscrew the 4 bolts (A).
3. Position the axle plate at the desired height.

**Note:** The thickest part of the axle plate must be at the bottom.


4. Retighten the bolts.
5. Then adjust the other axle plate.

## Settings and adjustments

**Note:** When changing the seat height, seat angle and axle plate position of the front wheel, the brake and the front wheel telescope must, if required, be readjusted.


### Setting the fork angle

If you do not raise/lower the axle plate and telescope proportionally, the fork angle changes with respect to the floor. To retain the favourable, manoeuvrable ride characteristics of the Roxx, the swivel castor axle must be perpendicular to the floor.

You should always set the swivel castor axles to be vertical with respect to the floor (this is 90°).

1. Loosen the socket screws (A) slightly so that the front wheel block (C) can be turned by hand.
2. Then place a spirit level against the vertical side of the front wheel block.
3. Now use a 24 mm spanner to turn the eccentric bolt (B) until the front wheel block is vertical.
4. Retighten the socket screws.

**Note:** As this is difficult to see, you can check it by turning the fork from front to back. If the swivel castor axle is vertical, there is no measurable height difference between the floor and the bottom of the telescopic tube.


## Settings and adjustments

Roxx in combination with Trippel telescope		
Required seat height	Tyre size	Position at seat angle 0°
38 cm	5" x 22"	T1-A / A1
39 cm	5" x 22"	T1-B / A2
40 cm	5" x 22"	T1-C / A2
	5" x 24"	T1-C / A1
	6" x 22"	T1-A / A2
	6" x 24"	T1-A / A1
41 cm	5" x 22"	T1-D / A3
	5" x 24"	T1-D / A2
	6" x 22"	T1-B / A3
	6" x 24"	T1-B / A2
42 cm	5" x 22"	T1-E / A3
	5" x 24"	T1-E / A2
	6" x 22"	T1-C / A3
	6" x 24"	T1-C / A2
43 cm	5" x 22"	T1-E / A3
	5" x 24"	T1-E / A2
	6" x 22"	T1-C / A3
	6" x 24"	T1-C / A2
44 cm	5" x 22"	T1-F / A4
	5" x 24"	T1-F / A3
	6" x 22"	T1-D / A4
	6" x 24"	T1-D / A3
45 cm	6" x 22"	T1-E / A4
	6" x 24"	T1-E / A3
46 cm	6" x 22"	T1-F / A5
	6" x 24"	T1-F / A4
	6" x 26"	T1-F / A3

## Settings and adjustments

**Adjusting seat height and seat angle**  
 You can adjust the seat height and if required the seat angle of your Roxx to various positions. To adjust the seat height, the axle plates can be adjusted to 6 heights.


### Positions of axle plate: A1 to A5


- A1: 1st and 2nd axle hole from the top
- A2: 2nd and 3rd axle hole from the top
- A3: 3rd and 4th axle hole from the top
- A4: 4th and 5th axle hole from the top
- A5: 5th and 6th axle hole from the top

### Hole positions for Trippel telescope (=T1): T1-A to T1-F


T1 = trippel telescope  
 T2 = standard telescope  
 T3 = standard telescope extended


T1-A: lowest position


T1-B: one but lowest position


T1-C: 1 hole visible

T1-D: 2 holes visible

T1-E: 3 holes visible


T1-F: 4 holes visible

## Settings and adjustments

Roxx in combination with Standard telescope (fork angle adjustable)		
Required seat height	Tyre size	Position for seat angle ~ 4°
44 cm	6" x 22"	A3 / T2-A
	6" x 24"	A2 / T2-A
45 cm	6" x 22"	A3 / T2-B
	6" x 24"	A2 / T2-B
46 cm	6" x 22"	A4 / T2-C
	6" x 24"	A3 / T2-C
47 cm	6" x 22"	A4 / T2-D
	6" x 24"	A3 / T2-D
48 cm	6" x 22"	A4 / T2-D
	6" x 24"	A3 / T2-D
49 cm	6" x 22"	A5 / T2-E
	8" x 22"	A5 / T2-A
	6" x 24"	A4 / T2-E
	8" x 24"	A4 / T2-A
50 cm	6" x 22"	A5 / T2-F
	8" x 22"	A5 / T2-B
	6" x 24"	A4 / T2-F
	8" x 24"	A4 / T2-B
	6" x 26"	A3 / T2-F
	8" x 26"	A3 / T2-B
51 cm	6" x 22"	A6 / T2-G
	8" x 22"	A6 / T2-C
	6" x 24"	A5 / T2-G
	8" x 24"	A5 / T2-C
	6" x 26"	A4 / T2-G
	8" x 26"	A4 / T2-C
52 cm	6" x 22"	A6 / T2-H
	8" x 22"	A6 / T2-D
	6" x 24"	A5 / T2-H
	8" x 24"	A5 / T2-D
	6" x 26"	A4 / T2-H
	8" x 26"	A4 / T2-D

T2-A: lowest position


T2-G: 5 holes visible

T2-B: one but lowest position


T2-H: 6 holes visible

T2-C: 1 hole visible


T2-I: 7 holes visible

T2-D: 2 holes visible

T3-A: 4 holes visible

T2-E: 3 holes visible

T3-B: 5 holes visible


T2-F: 4 holes visible

## Settings and adjustments

Required seat height	Tyre size	Position for seat angle ± 4°
52 cm	8" x 26"	A4 / T2-D
53 cm	6" x 22"	A6 / T2-H
	8" x 22"	A6 / T2-D
	6" x 24"	A5 / T2-H
	8" x 24"	A5 / T2-D
	6" x 26"	A4 / T2-H
54 cm	8" x 26"	A4 / T2-D
	6" x 24"	A6 / T3-A
	8" x 24"	A6 / T2-E
	6" x 26"	A5 / T3-A
55 cm	8" x 26"	A5 / T2-E
	6" x 24"	A6 / T3-B
	8" x 24"	A6 / T2-F
56 cm	6" x 26"	A5 / T3-B
	8" x 26"	A5 / T2-F
	6" x 24"	A6 / T3-B
56 cm	8" x 26"	A6 / T2-G
57 cm	8" x 26"	A6 / T2-H
58 cm	8" x 26"	A6 / T2-H
59 cm	8" x 26"	A6 / T3-A
60 cm	8" x 26"	A6 / T3-B

### Positions of axle plate: A1 to A6

- A1: 1st and 2nd axle hole from the top
- A2: 2nd and 3rd axle hole from the top
- A3: 3rd and 4th axle hole from the top
- A4: 4th and 5th axle hole from the top
- A5: 5th and 6th axle hole from the top
- A6: 6th and 7th axle hole from the top


## Settings and adjustments

### Ride characteristics

The ride characteristics of the Roxx are adjustable.

1. Remove the wheel.
2. Unscrew the nut (A) on the inside.
3. Move the axle bearing in the wheel plate forward or backward (the last axle hole provides more stability, the front axle hole provides more comfort).
4. Retighten the nut (A).
5. Then move the other axle bearing.
6. Readjust the brakes.


If the axle bearing has already been mounted in the front position and the user wants an even more active adjustment, then it is possible to exchange the left and right axle plate.

**Note:** If this is done, then the foldable frame can no longer be fully folded.

### Brakes

When changing the seat height, the seat angle and the ride characteristics, it is necessary to readjust the brakes.

1. Unscrew the socket screw (B) at the bottom of the brake clamp.
2. Slide the brake to the desired position.
3. Retighten the socket screw (B) at the bottom of the brake clamp.

### Backrest angle (once-only adjustable)


1. Unscrew the socket screws (A) in the backrest pivot points a number of turns.
2. Set the backrest to the required angle.
3. Then retighten the socket screws.

**Note:** Ensure that the left and the right push-bars are at the same angle using the graduation marks.

### Backrest height

The backrest can be set to different heights.

1. Remove the backrest.
2. Unscrew the socket screws (or wing nuts) (B).
3. Remove the backrest tubes.
4. Unscrew the socket screws (C) with the nut.
5. Select the desired backrest height.
6. Retighten the socket screws and nut.
7. Refit the backrest tubes on to the seat frame.
8. Retighten the socket screws (or wing nuts) (B).


**Note:** Ensure that the left and right backrest tube are at the same height.

**Note:** Check whether the backrest is securely fitted before sitting in the wheelchair.

### Adjusting the push-handle

1. Unscrew the wing nuts (A) one turn.
2. Position the push-handle at the desired height, the heights are indicated by depressions.
3. Retighten the wing nuts.

## Settings and adjustments

### Armrests


The armrests can be adjusted in height, depth and width.

#### Setting the height

1. Unscrew the socket screw (B) at the front of the armrest a couple of turns, or pull the knob.
2. Position the armrest at the required height.

**Note:** The armrest must click into position before the screw is retightened/ the knob is pressed.

3. Retighten the socket screw or press in the knob.

#### Adjusting the depth and width of the PUR armrest

1. Unscrew the screws (C).
2. Adjust to the desired width and/or depth.
3. Retighten the screws.

#### Adjusting the depth and width of the padded armrest


1. Unscrew the screws completely (different hole in armrest) or partially (using a spanner).
2. Adjust to the desired width and/or depth.
3. Retighten the screws

### Legrests


The legrests can be adjusted to the dimensions and wishes of the user.

#### Adjusting the lower leg length

1. Unscrew the socket screw (A) one turn.
2. Slide the legrest in or out to obtain the desired lower leg length. Use the graduation marks on the legrest to determine the exact length.


## Settings and adjustments

**Note:** Ensure that the legrest is at least 5 cm from the floor.

3. Retighten the socket screw.


#### Locking the position of the ankle hinge

1. Adjust the ankle angle to the required position.
2. Drill a hole in the hinge.
3. Lock the ankle hinge using the locking pin provided.

### Anti-tipping supports

The anti-tipping supports can be adjusted using the lever (A). The recommended height is 3 cm from the floor.

#### Please note:

- Make sure that the pin clicks properly into the opening.
- Inform the user if the anti-tip support is folded away.

#### Tip:

- The sliding sleeve can be used to mark a standard position.

The anti-tip support - if properly adjusted - protects the wheelchair user from tipping backwards. The risk of tipping the wheelchair backwards is influenced by:

- the position of the user's center of gravity in the wheelchair (which is determined by factors such as axle position, backrest angle, seat angle, seat tilt, and weight distribution over the body (e.g., amputations).
- any restlessness or spasms of the user.
- the slope on which you are driving.

Adjust the anti-tip as necessary based on the above considerations.


## Folding the wheelchair

**Folding the wheelchair** (only if the wheelchair has a foldable inner frame)

*Removing the backrest from the backrest frame*

1. Lift the left and right knobs (A) of the top fittings at the same time.
2. Pull the backrest forward so that the U profile releases from the backrest attachment points.
3. Release the knobs of the top fittings.
4. Lift out the backrest.


*Removing the seat*

1. Pull up the seat at the back.
2. Lift the seat out of the frame.

*Removing the legrests*

1. Lift the legrest out of the adapter, or loosen the star knob/ socket screw (B) that is used to lock the legrest in the seat frame and remove it by moving it forward.


*Standard side fitted with armrest*

1. Unlock the side.
2. Remove the side.

*Removing the headrest (if present)*

1. Loosen the star knob (C).
2. Remove the headrest from the adapter.


## Folding the wheelchair

*Backrest tubes*

1. Unscrew the socket screws (A) or wing nuts.
2. Remove the backrest tubes.

*Folding*

The chair is now ready to be folded.

1. Pull the middle round tube upwards.
2. The chair will fold.

 **Note:** Ensure that your fingers do not get caught between the parts of the frame.

 **Note:** Ensure that the chair is completely unfolded before using it.

 **Note:** Check that all parts are properly fixed before you sit in your wheelchair.

## Use

## Use

The Roxx is a wheelchair for users with various handicaps. Transfers from the wheelchair to, for example, a bed, toilet or shower and vice versa can be made in various ways. The extent to which the wheelchair occupant is able to cooperate determines how he or she should be lifted. Both active and passive transfers are possible with the Roxx, if required with the aid of the care provider or using a lift.

### Getting in and out of the chair yourself

 **Note:** Engage the brake on the wheelchair before getting in or out.

 **Note:** Do not use the brake as an aid to get up.

- Fold up the footrests.
- Lift the legrest up a little.
- Swing the legrest outward\*, or remove it.


If you have a Comfort legrest, fold the calfrest backwards if necessary.

- Lift your leg off the calfrest.
- Press the round knob so that the calfrest folds backwards.
- Take care when getting in and out of the chair.


In the event of a sideward transfer, you can, if required, remove or if possible fold up the sides.

### Tackling obstacles


When tackling obstacles on your own, Life & Mobility recommends practicing extensively under the supervision of a therapist.

*\* Important:  
take care not to disengage the  
band brake while folding away  
the leg support.*

Kerb

Gentle slope

 **Note:** When travelling up a slope, the centre of gravity will move backwards, making it easier for the wheelchair to tip over backwards.

Life & Mobility advises less experienced users to adjust the wheelchair somewhat so that its characteristics are less active and to fit it with an anti-tipping support.

Tackling obstacles with an accompanying person:

- Push down with your foot on the tipping lever and push down on the push-handles.
- Push the wheelchair onto the obstacle and remove your foot from the tipping lever.
- Push the wheelchair further and lift the back wheels over the obstacle using the push-handles.

You can mount a kerb forwards, but it is better to descend backwards.

 **Note:** Keep your speed under control on slopes.

### Brakes

The brakes of the Roxx are correctly adjusted as standard.

- The correct manner of braking is by pushing both the left and the right brake levers forward.
- The brake decouples by exerting a slight backward pressure on the brake lever. The brake will be returned to the correct position by a spring.
- The brakes are only intended as a parking brake, not to slow down the wheelchair.

## Transport


## Transport

### Taxibus

The Roxx was designed to be secured in a forward direction in a taxi van. The Roxx was tested positively as a crash-test safe wheelchair in accordance with ISO 7176-19:2008. It can therefore be transported with the occupant in it in a taxibus (see also the sticker on the wheelchair).

However, whether the wheelchair occupant can be safely transported depends on the state of his or her health; the user's physician should be consulted. Life & Mobility recommends transferring the occupant from the Roxx into a fixed seat if possible.

- The Roxx has been tested, with the test dummy fixed in position using a 3-point safety belt. The suitability of the Roxx in crashes or other types of belts other than those specified in the above standard has not been tested.
- The Roxx was tested in various configurations.
- Life & Mobility recommends that the chair should always be fitted with a Matrixx headrest.
- For preventive safety reasons, Life & Mobility installs, ex factory, a headrest adapter on the backrest.
- The Roxx must be fastened to the floor of the taxibus using an approved\* four-point securing system.

**Note:** If there is no ex factory crashtest safe sticker on the wheelchair, it is not a wheelchair that can be transported safely with the occupant in it.

**Note:** Only lift the wheelchair using fixed frame components and not removable accessories such as armrests, legrests seating, back support etc.

**Note:** The Roxx is only crashtest safe in combination with the original Matrixx elements and parts.

\* In accordance with EN10542


Figure 1


Figure 2


Figure 3


Figure 4

- The positions of the four securing points on the Roxx are shown with a hook symbol sticker (see Figure 1).
- Easy access and manoeuvrability in motor vehicles can be significantly influenced by wheelchair size and agility. Smaller wheelchairs with a shorter turning circle generally offer better accessibility to the vehicle and manoeuvrability to a forward position. The Roxx was assessed positively for securability according to ISO 7176-19:2008, Annex D.
- The Roxx must be transported with the occupant facing forward. The wheelchair must be secured in accordance with the instructions of the manufacturer of the securing system.
- The Roxx should be fastened using anchoring systems that are suitable for the 4 attachment points).
- The wheelchair occupant must wear an approved\* three-point safety belt in the taxibus. These safety belts must be attached to the taxibus.
- The lap belt must be positioned low across the occupant's waist. The angle between the belt and the floor must be between 30° and 75° (a larger angle is better) (see Figure 2).
- The lap belt must run over the front of the body following the line where the hip and the pelvis meet. The belt must closely follow the body here. The top belt must fully support the occupant as illustrated in Figure 3.
- The diagonal belt must be positioned across the occupant's shoulder and chest as shown in Figure 4.
- The belt must be sufficiently tight but should not cause the occupant discomfort.
- The belt must not be twisted.

## Transport

- The belt must not come into contact with components of the wheelchair (such as armrests or wheels) in such a way that it does not fully support the occupant's body (see Figure 5).
- The approved three-point safety belt is intended to prevent injury to the head and chest.
- If a worktop is fitted, it must be removed in order to prevent injury to other passengers. The worktop should be safely stowed away elsewhere in the vehicle.
- If possible, any accessories must be secured or removed.
- Following an accident, a representative of the factory must inspect the wheelchair to determine whether it can still be used.
- The attachment points, structural components and frame components must not be modified or replaced without first consulting the manufacturer.

### Seat belt

The seat belt is intended to keep the occupant in position. Attach it to or around the metal seat frame.

### Car

In addition to a taxibus, the Roxx can also be transported folded up in the boot of a passenger car. See page 22 and 23 for instructions on how to fold up the wheelchair.

Figure 5

## Transport

### Warnings

- In accordance with the agreement, the Roxx base was declared collision-test safe (ISO 7176-19:2008) in combination with specific orthoses and was issued with CE marking.
- The Roxx has been dynamically tested in a forward direction, where the test dummy (= 103 kg) was secured with a hip belt and a shoulder belt.
- Secure the chair in a taxi van with a 4point fastening system. Also use a 3-point safety belt for transporting the wheelchair and occupant to reduce the risk of head and chest injuries.
- The sitting posture belt is unsuitable as a safety belt.
- Worktops should be removed and stored safely in the taxi.
- All accessories should be secured or removed and safely stored in the taxi.
- After every accident, the manufacturer's representative should examine the wheelchair to determine whether it can be used again.
- The fastening points or structural and frame components should never be changed and/or replaced without first consulting the manufacturer.
- Any batteries should be "spill-proof".
- The belt clasp must be positioned so that no wheelchair parts can come in contact with it during a collision.

## Maintenance

### Maintenance

*Minor maintenance (every two weeks)*

You can do minor maintenance on the Roxx yourself. Clean your wheelchair with a damp cloth, water and an all-purpose cleaner.

**Note:** Never use abrasive and aggressive cleaning agents or solvents including thinner, dry cleaning fluid or similar.

**Note:** Over-tightening the screws can damage the wheelchair.

**Tip:** Use Matrixx Cleaner for cleaning the upholstery.

### *Tyres and wheels*

Check the tyre pressure and condition of the tyres (for ageing) regularly. To propel the wheelchair comfortably, meaning with the lowest possible rolling and turning resistance, it is important for the tyres to be at the correct pressure. The wheelchair will be more difficult to propel and manoeuvre if the tyre pressure is too low. Use, if necessary, a valve nipple for pumping up the tyres (for example, when using a bicycle pump). The maximum tyre pressure is shown on the tyres. Life & Mobility recommends a maximum tyre pressure of 6 bar.

**Note:** Tyre pressures that are too low will have a negative impact on the braking power.

Check the wheel tension on a regular basis (once per quarter, depending on the intensity of use) to make sure that your wheelchair continues to ride comfortably and safely. Ask your dealer for assistance if necessary.

## Maintenance

### Major maintenance (annual)

Your dealer will be happy to carry out the major maintenance of your wheelchair. You can then be sure that you will be getting the best possible use out of it. You should also consult your dealer if components need replacing.

## Technical problems

### Technical problems

In the unlikely event that you experience technical problems with your wheelchair, please contact your dealer. The dealer will decide whether he can repair the chair or whether it will have to be returned to the factory. You may be able to solve some problems yourself.

*Parts that can be set or adjusted, including push-handles, armrests, legrests, etc. cannot be moved*

- These components can become jammed in the frame if too much pressure is exerted on them. Carefully tap the components loose with a plastic hammer.

*Socket screws used to adjust the wheelchair cannot be loosened*

- The screws may have been over-tightened. Try placing a tube over the long end of the Allen key to extend it. Now try loosening the screw.

**Note:** Do not do this to tighten a screw. The screw can be over-tightened and cause damage.

- The socket head may be stripped: Try gripping the screw with a monkey wrench. If this fails consult your dealer.

*The brakes have insufficient braking power or do not brake with similar braking power on both sides*

- You should readjust the brakes, see the description on page 18.

Regularly check the tyre pressure, and if required, pump up the tyres.

## Guarantee

### Terms and Conditions of Warranty

The Life & Mobility products are made with care and are carefully inspected before leaving the factory. Should it emerge that a product does not meet expectations please contact the party where you purchased the product.

Life & Mobility respects the laws and regulations that protect the consumer and that apply in the country in which the product is purchased.

The warranty does not apply to defects caused by: failure to observe the instructions for use and maintenance, use other than normal use, wear and tear, negligence, overloading, third-party accidents, non-original parts used and defects not caused by the product.

This warranty replaces any other warranty, either stipulated by law or communicated verbally, except for anything guaranteed in writing by Life & Mobility.

Warranties are only valid within the EU.

## Used wheelchairs/environment

### Disposal of the product and packaging

From the environmental viewpoint, we recommend disposing of the packaging and the wheelchair in such a way that the materials used can be efficiently recycled.

In the majority of cases, you can hand in your used wheelchair at your wheelchair dealer. If this is impossible, you can hand in the wheelchair at your local refuse collection location.

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