

PRESTO



Out and about
with confidence


Inspired by life




PRESTO

Out and about with confidence

If you are looking for a safe and first-class mobility scooter with excellent driving comfort thanks to the unique INDEGO suspension system, then the Presto is the model for you. Its tight turning circle makes it ideal for both outdoor and indoor use. When it comes to the design, we have carefully considered all the safety aspects associated with using a mobility scooter. For example, the LED lights ensure you are highly visible in traffic. The clear and well thought-out control panel means you can set off effortlessly. A quick trip to the shops or just going out to see family or friends. Its beautiful and well thought-out design means you can get out and about with confidence!



15 km/h



44 km

INDEGO
suspension system



SDS SAFE DRIVING SOLUTIONS

Safety

- To increase your visibility, the indicators are positioned high up and the rear lights are designed so they can also be seen from the side.
- Directional headlight with LED bulbs for better vision when cornering.
- You can keep your eyes on the road because the control buttons can be identified by touch and are positioned within thumb's reach.
- When cornering, your speed is automatically restricted to a safe level so you remain in full control.
- The indicators switch off automatically after a turn and can be operated safely without having to let go of the steering wheel.

EDS ERGONOMIC DRIVING SOLUTIONS

Comfortable ride

- Outstanding ride comfort thanks to the unique INDEGO suspension system.
- Clearly legible pictograms and information thanks to a large and well laid-out control panel.
- Optimum sitting position thanks to the multi-adjustable comfortable seat: height, depth, backrest angle, armrest height, width and seat cushion angle are fully adjustable.
- Generous leg and foot room thanks to the ergonomic design.

Features

- Adjustable speed up to 15 km/h.
- Maximum operating range of 44 km*.
- Tight turning circle of 160 cm.
- High load capacity of 160 kg.
- Removable plastic shopping basket.
- Innovative tubeless tyres with low rolling resistance.

* Under ideal circumstances



Digital display



Plenty of foot room



Directional headlight



LED lights