


Vivo



The little scooter that  
expands your world!





# Vivo

## Demountable travel scooter with excellent suspension

Embrace the freedom offered by the Vivo. The Vivo travel scooter is easy to disassemble so you can take it with you wherever you go. Its very small turning radius makes the Vivo easy to manoeuvre, while its wide tyres give it an extremely high degree of stability. Its front and rear suspension provide unique passenger comfort and, with an operating range of 28 km, you can leave home without a worry. Vivo, the little scooter that expands your world!

### Specifications

- **Comfortable thanks to wide tyres and suspension on the rear axle**
- **Highly manoeuvrable due to very small turning radius**
- **Large operating range; up to 28 km\***
- Max. carrying weight of 120 kg
- Plenty of foot space
- Batteries (22 Ah) and charger (3 A) included

\* Under ideal conditions


### Adjustability

- **Can be dismantled in at least 5 components for easy transportation and storage**
- Adjustable and retractable handlebars
- Comfortable, removable swivel seat with adjustable seat height
- Foldable armrests with adjustable angle
- Removable shopping basket
- Removable battery set that can be separately charged
- Freewheel lever
- **Easy colour personalisation by change of shock-resistant bonnet set**



### Display

- Simple control panel with magnetic lock, throttle levers, horn, battery indicator and speed adjustment
- Adjustable speed (max. 8 km/h)



### Safety

- Anti-tip wheels included

### Available colours



### Options:

- Crutch holder
- Rear-view mirror
- Backrest bag
- Seatbelt
- Cover